

Therabody

RecoveryAir PRO

Compression Half Jacket

Table of Contents

EN	1-2
DE	3-4
FR	5-6
FRCA	7-8
IT	9-10
ES	11-12
NL	13-14
PL	15-16
NO	17-18
SV	19-20
RU	21-22
CN	23-24

RecoveryAir PRO Compression Half Jacket

Instructions For Use

Before first use, read the [RecoveryAir PRO User Manual](#) for instructions on operating the full pneumatic compression system and for the full Safety Warnings and Precautions.

Wearing Instructions

The RecoveryAir PRO Compression Half Jacket can be used lying down, or in a seated position. The RecoveryAir PRO Compression Half Jacket is uniquely designed to treat the arms, shoulders, back, chest, and core. The arm can be treated held straight out to the side (90 degrees), down in a natural position, or held supported using the optional strap.

1. Place the RecoveryAir PRO pump where you can easily reach it during treatment.
2. Attach the garment hose to one of the air outlets on the pump and cap the other air outlet with the plug provided with the RecoveryAir PRO pump.
3. Lay RecoveryAir PRO Compression Half Jacket on a flat surface.
4. Zip the fastener to the vest, with the inner lining facing up (see Figure 1).
5. Open the zipper at the wrist of the arm sleeve. This allows you to use both hands to put on the garment.
6. Loosen the side straps on the fastener until you nearly reach its full length. This allows you to grasp them after you put the garment on.
7. Put on the garment, as you would a jacket. Make sure the armpit chamber is tucked into the arm sleeve.
8. Pull the side fastener around your torso. Use the zipper attached to the straps to connect with the main section of the garment. Zip the wrist zipper closed. Note: Make sure the padding is not folded and fits comfortably around the torso.
9. Tighten the straps at the side for a comfortable fit. If you experience discomfort during treatment, adjust the buckles at the side to loosen the garment.


Figure 1 – Assembling the RA Half Jacket

Using the Optional Arm Support Strap

1. Before turning on the RecoveryAir PRO pump power button, fasten the strap to the torso.
2. Then, extend the strap to its full length by pulling up on the buckle (see Figure 2) and fasten the strap to the arm sleeve.
3. Hold your arm in a comfortable position, but not bent, to allow the garment to fill with air and then operate the pump for one inflation cycle.
4. When the pump deflates, bend your arm into the desired position.
5. Adjust the strap to the desired tension by pulling on it with your free hand.

Operating the RecoveryAir PRO Pump

Select a pressure of up to 50 mmHg for use with RecoveryAir PRO Compression Half Jacket.

Note: Filling the garment may take longer on the first inflation. The garment should not fit too tightly. You should be able to insert three fingers into the garment when it is fully closed before inflation.

Caution

- Do not stand or walk while wearing the garment. Put on the garment while in a seated position and remove it before rising from the treatment position.
- Do not use the RecoveryAir PRO Compression Half Jacket at pressures over 50 mmHg. Using the RecoveryAir PRO Compression Half Jacket at pressure above 50 mmHg may cause discomfort and may reduce the effect of treatment.
- Do not inflate the RecoveryAir PRO Compression Half Jacket without wearing it over the intended body area ("on empty") or with open zippers. Doing so can damage the garment.
- Do not apply excessive force to the RecoveryAir PRO Compression Half Jacket straps. Do not use the straps for any purpose other than intended by Therabody. Use of excessive force and/or misuse shall void the warranty.
- Stop treatment and consult your licensed medical practitioner if you experience the onset of pain, bruising or irritation during or after treatment.
- Do not use the RecoveryAir PRO Compression Half Jacket if you are pregnant.

Warranty

For full warranty information, please visit www.therabody.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address:

Therabody - Warranty
Attn: Customer Service,
6100 Wilshire Blvd., Suite 200
Los Angeles, CA 90048

Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.


Figure 2 - Extending the Optional Arm Strap

RecoveryAir PRO Compression Half Jacket

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