
Therabody

RecoveryAir PRO

Compression Pants

Table of Contents

EN	1-2
DE	3-4
FR	5-6
FRCA	7-8
IT	9-10
ES	11-12
NL	13-14
PL	15-16
NO	17-18
SV	19-20
RU	21-22
CN	23-24

RecoveryAir PRO Compression Pants

Instructions For Use

Before first use, read the [RecoveryAir PRO User Manual](#) for instructions on operating the full pneumatic compression system and for the full Safety Warnings and Precautions.

Wearing Instructions

Use the RecoveryAir PRO Compression Pants in a supine position (lying face up on a bed or couch). Always wear light, comfortable clothing under the RecoveryAir PRO Compression Pants (sweatpants or leggings and socks, free of buttons or zippers that may rub against the skin) to absorb perspiration, prevent irritation, and keep your RecoveryAir PRO Compression Pants clean. Do not wear your RecoveryAir PRO Compression Pants over bare skin.

1. Position the RecoveryAir PRO pump so that you can reach it easily from the massage position. Lay the RecoveryAir PRO Compression Pants on the bed or couch, with the logo side facing up.
2. Open the zippers on both sides of your RecoveryAir PRO Compression Pants halfway. Get into the Pants and close the zippers all the way to the top.
3. To increase contact with your abdomen, you may place a pillow or bolster under your knees.

Note: The RecoveryAir PRO Compression Pants should not fit tightly. A pair of Expanders are available for an increase in width.

Using the Optional Foot Straps

You may use the Foot Straps to form a customized “boot” for your foot (see Figure 1).

1. Adjust the lengthwise straps at the tension buckle. Located at the low section of the legs, pull the open side of the strap to raise the boot to an angle that is comfortable. To lower the boot angle back to parallel, loosen the strap at the buckle.
2. After a few inflation cycles you may adjust the angle of the boots by releasing the buckle or pulling on the strap.
3. When the massage is complete, leave the straps in place for your next massage.


Figure 1

Using the Optional Waist and Leg Extenders

You may use the Optional Expanders to increase the width or circumference of the waist and legs. (see Figure 2).

1. Completely open the zippers located at the sides of the pants.
2. Then, zip the respective expanders (1 for each side, Right or Left) all the way from the feet to the waist/torso.
3. Once you feel the pants are adjusted to a comfortable fit, you can begin your massage.

Operating the RecoveryAir PRO Pump

1. Securely insert the garment hose bundle connectors, with the logo facing up, into the air sockets on the pump.
2. Select a pressure of up to 50 mmHg for use with your RecoveryAir PRO Compression Pants.

Note: Filling the garment may take longer on the first inflation. The garment should not fit too tightly. You should be able to insert three fingers into the garment when it is fully closed before inflation.

Caution

- Do not use the RecoveryAir PRO Compression Pants at pressures over 50 mmHg. Using the RecoveryAir PRO Compression Pants at pressure above 50 mmHg may cause discomfort and may reduce the effect of treatment.
- Do not stand or walk while wearing the garment. Put on the garment while in a seated position and remove it before rising from the treatment position.
- RecoveryAir PRO Compression Pants are not intended for use in a standing or seated position. Use in a reclined or supine (lying face up) position only.
- Do not inflate the pants without wearing it over the intended body area ("on empty") or with open zippers. Doing so can damage the garment.
- Do not apply excessive force to the pant straps. Do not use the straps for any purpose other than intended by Therabody. Use of excessive force and/or misuse shall void the warranty.
- Stop treatment and consult your licensed medical practitioner if you experience the onset of pain, bruising or irritation during or after treatment.
- Do not use the RecoveryAir PRO Compression Pants if you are pregnant.

Warranty

For full warranty information, please visit www.therabody.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address:

Therabody - Warranty
Attn: Customer Service,
6100 Wilshire Blvd., Suite 200
Los Angeles, CA 90048

Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.


Figure 2

RecoveryAir PRO Compression Pants

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