
Therabody

RecoveryAir

Compression Sleeve

Table of Contents

EN	1
DE	2
FR	3
FRCA	4
IT	5
ES	6
NL	7
PL	8
NO	9
SV	10
RU	11
ZHCN	12

RecoveryAir Compression Sleeve

Instructions For Use

Before first use, read the RecoveryAir User Manual for instructions on operating the full pneumatic compression system and for the full Safety Warnings and Precautions.

Wearing Instructions

The RecoveryAir Compression Sleeve can be used lying down, or in a seated position. The RecoveryAir Compression Sleeve is uniquely designed to treat the full arm, from hand to shoulder. The arm can be treated held straight out to the side (90 degrees) or down in a natural position. The RecoveryAir Compression Sleeve is designed to fit on either the right or the left arm.

1. Insert your arm into the RecoveryAir Compression Sleeve
2. Pull the RecoveryAir Compression Sleeve over your shoulder, fully covering the shoulder.
3. Then, wrap the fastening strap around your body, under the untreated arm.
4. Fasten the strap and adjust the buckle for a comfortable fit. If you experience discomfort during treatment, adjust the buckles at the side to loosen the garment.

Note: Filling the garment may take longer on the first inflation. The garment should not fit too tightly. You should be able to insert three fingers into the garment when it is fully closed before inflation.

Caution

- Do not stand or walk while wearing the RecoveryAir Compression Sleeve. Put on the garment while in a seated position and remove it before rising from the treatment position.
- Do not inflate the RecoveryAir Compression Sleeve without wearing it over the intended body area ("on empty"). Doing so can damage the garment.
- Do not apply excessive force to the RecoveryAir Compression Sleeve straps. Do not use the straps for any purpose other than intended by Therabody. Use of excessive force and/or misuse shall void the warranty.
- Stop treatment and consult your licensed medical practitioner if you experience the onset of pain, bruising or irritation during or after treatment.
- Do not use the RecoveryAir Compression Sleeve if you are pregnant.

Warranty

For full warranty information, please visit www.therabody.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address:

Therabody - Warranty
Attn: Customer Service
6100 Wilshire Blvd., Suite 200
Los Angeles, CA 90048

Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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